

## Product Datasheet

# Recombinant Human TNFRSF10C C-fc Protein, Fc & His Tag orb1551515

|                        |  |
|------------------------|--|
| <b>Description</b>     | Recombinant Human TNFRSF10C C-fc Protein, Fc & His Tag   |
| <b>Species/Host</b>    | Mammal   |
| <b>Reactivity</b>      | Human  |
| <b>Endotoxins</b>      |  |
| <b>Target</b>          | TNFRSF10C  |
| <b>Note</b>            | For research use only.   |
| <b>Purity</b>          | ≥90% as determined by SDS-PAGE   |
| <b>Source</b>          | This protein contains the human TNFRSF10C(Ala26-Ala221) was fused with the C-terminal Fc Tag and expressed in Mammalian cells. |
| <b>Uniprot ID</b>      | <a href="#">O14798</a>   |
| <b>Expiration Date</b> | 6 months from date of receipt.   |