

Product Datasheet

Recombinant Human TNFRSF11B Protein, His Tag orb1551513

Description	Recombinant Human TNFRSF11B Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	
Target	TNFRSF11B
Note	For research use only.
Purity	?90% as determined by SDS-PAGE
Source	This protein contains the human TNFRSF11B(Glu22-Leu401) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	O00300
Expiration Date	6 months from date of receipt.