

Product Datasheet

Recombinant Human TNFRSF18/GITR Protein, Fc Tag orb1551510

Description	Recombinant Human TNFRSF18/GITR Protein, Fc Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	
Target	TNFRSF18
Note	For research use only.
Purity	?90% as determined by SDS-PAGE
Source	This protein contains the human TNFRSF18(Gln26-Gln161) was fused with the C-terminal Fc Tag and expressed in Mammalian cells.
Uniprot ID	Q9Y5U5
Expiration Date	6 months from date of receipt.