

Product Datasheet

Recombinant Mouse Tnfrsf1a Protein, No Tag orb1551508

Description	Recombinant Mouse Tnfrsf1a Protein, No Tag
Species/Host	E. coli
Reactivity	Mouse
Endotoxins	Please contact our customer service for more information.
Target	Tnfrsf1a
Note	For research use only.
Purity	?90% as determined by SDS-PAGE
Source	This protein contains the mouse Tnfrsf1a(Ile22-Ala212) was fused without Tag and expressed in E. coli.
Uniprot ID	P25118
Expiration Date	6 months from date of receipt.