

Product Datasheet

Recombinant Human TNFRSF1B Protein, His Tag orb1551507

Description	Recombinant Human TNFRSF1B Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	
Target	TNFRSF1B
Note	For research use only.
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human TNFRSF1B(Lys288-Ser461) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	P20333
Expiration Date	6 months from date of receipt.