

Product Datasheet

POU5F1 Antibody orb1562902

Description	POU5F1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ChIP, ICC, IHC-P, IP, WB
Immunogen	A synthesized peptide derived from human Oct4
Target	POU5F1
Form/Appearance	Liquid
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Monoclonal
MW	Calculated MW: 39 kDa; Observed MW: 45 kDa
Clone ID	R08-6K5
Uniprot ID	Q01860
NCBI	5460
Dilution Range	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/50 ChIP: 1/20
Expiration Date	12 months from date of receipt.