

Product Datasheet

FEN1 antibody orb49073

Description	Rabbit polyclonal antibody to FEN1
Reactivity	Human, Mouse, Rat
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	Recombinant fusion protein of human FEN-1(NP_004102.1).
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	48kDa
Uniprot ID	P39748
NCBI	2237
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:10-100
Expiration Date	12 months from date of receipt.