

Product Datasheet

Human SAA4 ELISA Kit [Out of stock] orb407360

Description	Human Serum amyloid A-4 protein(SAA4) ELISA kit
Reactivity	Human
Range	15.6 ng/ml-1000 ng/ml
Sensitivity	3.9 ng/ml
Target	SAA4
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P35542
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm