

Product Datasheet

Human SAA1 ELISA Kit [Out of stock] orb407359

Description	Human Serum amyloid A,SAA ELISA Kit
Reactivity	Human
Range	156 ng/mL-10000 ng/mL
Sensitivity	39 ng/mL
Target	serum amyloid A1
Note	For research use only.
Sample Types	serum, urine
Uniprot ID	P0DJ18
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, urine• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm