

Product Datasheet

Human BSA ELISA Kit [Out of stock] orb406187

Description	Human bovine serum albumin (BSA) ELISA kit
Reactivity	Human
Range	0.312 ?g/mL-80 ?g/mL
Sensitivity	0.312 ?g/mL
Target	bovine serum albumin (BSA)
Note	For research use only.
Sample Types	serum, plasma
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma• Assay Time: 1-5 hours• Sample Volume: 50-100 ?l• Detection Wavelength: 450 nm