

Product Datasheet

Human SAA2 ELISA Kit [Out of stock] orb405314

Description	Human serum amyloid A2 (SAA2) ELISA kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Sensitivity	0.078 ng/mL
Target	serum amyloid A2 (SAA2)
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates. • Assay Time: 1-5 hours • Sample Volume: 50-100 µL • Detection Wavelength: 450 nm