

Product Datasheet

LpxT Peptide orb1787771

Description	LpxT Peptide
Tested Applications	Blocking
Target	lpxT
Form/Appearance	Liquid
Concentration	200 ug/mL
Preservatives	PBS pH 7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide
Storage	Store LpxT peptide at -20°C, stable for one year.
Note	For research use only.
Uniprot ID	P76445
NCBI	946693 , 2 , 76445
Expiration Date	6 months from date of receipt.
Application Notes	The peptide is used for blocking the activity of the LpxT antibody.