

Product Datasheet

Human PDCD6IP ELISA Kit orb407231

Description	Human Programmed cell death 6-interacting protein(PDCD6IP) ELISA kit
Reactivity	Human
Range	47 pg/mL-3000 pg/mL
Sensitivity	11.7 pg/mL
Target	PDCD6IP
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q8WUM4
Expiration Date	Please enquire.
Application Notes	Detection Wavelength: 450 nm