

## Product Datasheet

### Human MAP1LC3B ELISA Kit orb407206

|                   |                                                                                   |
|-------------------|-----------------------------------------------------------------------------------|
| Description       | Human Microtubule-associated proteins 1A/1B light chain 3B(MAP1LC3B)<br>ELISA kit |
| Reactivity        | Human                                                                             |
| Range             | 18.75 pg/mL-1200 pg/mL                                                            |
| Sensitivity       | 4.67 pg/mL                                                                        |
| Target            | MAP1LC3B                                                                          |
| Note              | For research use only.                                                            |
| Sample Types      | serum, plasma, tissue homogenates, cell lysates                                   |
| Uniprot ID        | <a href="#">Q9GZQ8</a>                                                            |
| Expiration Date   | Please enquire.                                                                   |
| Application Notes | Detection Wavelength: 450 nm                                                      |