

## Product Datasheet

### Human CD55 ELISA Kit orb407090

<b>Description</b>	Human decay-accelerating factor(DAF/CD55)ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.39 ng/ml- 25ng/ml
<b>Sensitivity</b>	0.09 ng/ml
<b>Target</b>	CD55
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">P08174</a>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	Detection Wavelength: 450 nm