

Product Datasheet

Human GADD45A ELISA Kit orb407075

Description	Human Growth arrest and DNA damage-inducible protein GADD45 alpha(GADD45A) ELISA kit
Reactivity	Human
Range	25 pg/mL-1600 pg/mL
Sensitivity	6.25 pg/mL
Target	GADD45A
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P24522
Expiration Date	Please enquire.
Application Notes	Detection Wavelength: 450 nm