

## Product Datasheet

### Human HSPA8 ELISA Kit orb406978

<b>Description</b>	Human Heat shock cognate 71 kDa protein(HSPA8) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	31.25 pg/mL-2000 pg/mL
<b>Sensitivity</b>	7.8 pg/mL
<b>Target</b>	HSPA8
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">P11142</a>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	Detection Wavelength: 450 nm