

## Product Datasheet

### Human HSP70 ELISA Kit orb406971

<b>Description</b>	Human Heat Shock Protein 70,HSP-70 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	1.56 ng/mL-100 ng/mL
<b>Sensitivity</b>	0.39 ng/mL
<b>Target</b>	HSPA1B
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">P0DMV9</a>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	Detection Wavelength: 450 nm