

## Product Datasheet

### Human NGB ELISA Kit orb406855

|                   |                                   |
|-------------------|-----------------------------------|
| Description       | Human Neuroglobin, NGB ELISA Kit  |
| Reactivity        | Human                             |
| Range             | 0.27 ng/ml - 200 ng/ml            |
| Sensitivity       | 0.27 ng/ml                        |
| Target            | NGB                               |
| Note              | For research use only.            |
| Sample Types      | serum, plasma, tissue homogenates |
| Uniprot ID        | <a href="#">Q9NPG2</a>            |
| Expiration Date   | Please enquire.                   |
| Application Notes | Detection Wavelength: 450 nm      |