

Product Datasheet

Human GAPDH/G3PDH ELISA Kit orb406296

Description	Human glyceraldehyde-3-phosphate dehydrogenase (GAPDH/G3PDH) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Sensitivity	7.81 pg/mL
Target	GAPDH/G3PDH
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P04406
Expiration Date	Please enquire.
Application Notes	Detection Wavelength: 450 nm