

Product Datasheet

Human SFRP5 ELISA Kit orb406261

Description	Human secreted frizzled-related protein 5(SFRP5) ELISA Kit
Reactivity	Human
Range	1.56 pg/mL-100 pg/mL
Sensitivity	0.39 pg/mL
Target	SFRP5
Note	For research use only.
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q5T4F7
Expiration Date	Please enquire.
Application Notes	Detection Wavelength: 450 nm