

## Product Datasheet

### Human PD1 ELISA Kit orb406125

|                   |                                                 |
|-------------------|-------------------------------------------------|
| Description       | Human Programmed Death 1(PD-1)ELISA KIT         |
| Reactivity        | Human                                           |
| Range             | 15.6 pg/mL-1000 pg/mL                           |
| Sensitivity       | 3.9 pg/mL                                       |
| Target            | PDCD1                                           |
| Note              | For research use only.                          |
| Sample Types      | serum, plasma, tissue homogenates, cell lysates |
| Uniprot ID        | <a href="#">Q15116</a>                          |
| Expiration Date   | Please enquire.                                 |
| Application Notes | Detection Wavelength: 450 nm                    |