

Product Datasheet

Human IL1 alpha ELISA Kit orb405978

Description	Human Interleukin 1?,IL-1? ELISA Kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Sensitivity	7.8 pg/mL
Target	IL1A
Note	For research use only.
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
Uniprot ID	P01583
Expiration Date	Please enquire.
Application Notes	Detection Wavelength: 450 nm