

## Product Datasheet

### Human Orexin A ELISA Kit orb405951

<b>Description</b>	Human Orexin A ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156 ng/mL-10 ng/mL
<b>Sensitivity</b>	0.039 ng/mL
<b>Target</b>	Orexin A
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">O43612</a>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	Detection Wavelength: 450 nm