

Product Datasheet

Human CTNNB1 ELISA Kit orb405928

Description	Human Catenin beta-1 (CTNNB1/CTNNB/OK/SW-cl.35/PRO2286) ELISA kit
Reactivity	Human
Range	15.6 pg/mL-1000 pg/mL
Sensitivity	3.9 pg/mL
Target	CTNNB1
Note	For research use only.
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
Uniprot ID	P35222
Expiration Date	Please enquire.
Application Notes	Detection Wavelength: 450 nm