

## Product Datasheet

### Human TLR4 ELISA Kit orb405856

|                          |   |
|--------------------------|---|
| <b>Description</b>       | human toll-like receptor 4,TLR4 ELISA Kit |
| <b>Reactivity</b>        | Human                                     |
| <b>Range</b>             | 0.156 ng/mL-10 ng/mL                      |
| <b>Sensitivity</b>       | 0.039 ng/mL                               |
| <b>Target</b>            | TLR4                                      |
| <b>Note</b>              | For research use only.                    |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates         |
| <b>Uniprot ID</b>        | <a href="#">O00206</a>                    |
| <b>Expiration Date</b>   | Please enquire.                           |
| <b>Application Notes</b> | Detection Wavelength: 450 nm              |