

Product Datasheet

Human p53 ELISA Kit orb405792

Description	Human p53/tumor protein,p53/TP53 ELISA Kit
Reactivity	Human
Range	9.38 pg/mL-600 pg/mL
Sensitivity	2.34 pg/mL
Target	TP53
Note	For research use only.
Sample Types	serum, cell culture supernates, urine, cerebrospinal fluid (CSF), tissue homogenates, cell lysates
Uniprot ID	P04637
Expiration Date	Please enquire.
Application Notes	Detection Wavelength: 450 nm