

Product Datasheet

Human ATG7 ELISA Kit orb405636

Description	Human Ubiquitin-like modifier-activating enzyme ATG7(ATG7) ELISA kit
Reactivity	Human
Range	28 pg/mL-1800 pg/mL
Sensitivity	7 pg/mL
Target	ATG7
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	O95352
Expiration Date	Please enquire.
Application Notes	Detection Wavelength: 450 nm