

Product Datasheet

Human VSIG4 ELISA Kit orb405547

Description	Human V-set and immunoglobulin domain-containing protein 4(VSIG4) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Sensitivity	7.81 pg/mL
Target	VSIG4
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9Y279
Expiration Date	Please enquire.
Application Notes	Detection Wavelength: 450 nm