

Product Datasheet

Human TNFAIP3 ELISA Kit orb405479

Description	Human Tumor necrosis factor alpha-induced protein 3(TNFAIP3) ELISA kit
Reactivity	Human
Range	23.5 pg/mL-1500 pg/mL
Sensitivity	5.8 pg/mL
Target	TNFAIP3
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P21580
Expiration Date	Please enquire.
Application Notes	Detection Wavelength: 450 nm