

Product Datasheet

Human TREM2 ELISA Kit orb405459

Description	Human Triggering receptor expressed on myeloid cells 2(TREM2) ELISA kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Sensitivity	0.078 ng/mL
Target	TREM2
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9NZC2
Expiration Date	Please enquire.
Application Notes	Detection Wavelength: 450 nm