

Product Datasheet

Human PCSK9 ELISA Kit orb405078

Description	Human Proprotein convertase subtilisin/kexin type 9(PCSK9) ELISA kit
Reactivity	Human
Range	0.45 ng/mL-30 ng/mL
Sensitivity	0.11 ng/mL
Target	PCSK9
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q8NBP7
Expiration Date	Please enquire.
Application Notes	Detection Wavelength: 450 nm