

Product Datasheet

Human WT1 ELISA Kit orb405052

Description	Human PRKC apoptosis WT1 regulator protein(PAWR) ELISA kit
Reactivity	Human
Range	25 pg/mL-1600 pg/mL
Sensitivity	6.25 pg/mL
Target	PAWR
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q96IZ0
Expiration Date	Please enquire.
Application Notes	Detection Wavelength: 450 nm