

Product Datasheet

P140 orb1821034

Description	HSPA8/HSC70 chaperone binder
Form/Appearance	Freeze dried solid
Storage	Store dry, frozen and in the dark
Note	For research use only.
Protein Sequence	RIHMOVYSKRS(PO3)GKPRGYAFIEY
Purity	> 95% by hplc
Formula	C ₁₁₁₇ H ₁₈₁ N ₃₂ O ₃₂ PS
MW	2638.94 Da
CAS Number	497156-60-2
Solubility (25°C)	Soluble in water
Expiration Date	6 months from date of receipt.