

## Product Datasheet

### Human ANGPT4 ELISA Kit orb407339

|                   |                                      |
|-------------------|--------------------------------------|
| Description       | Human Angiopoietin 4,ANG-4 ELISA Kit |
| Reactivity        | Human                                |
| Range             | Request Information                  |
| Sensitivity       | Request Information                  |
| Target            | ANGPT4                               |
| Note              | For research use only.               |
| Uniprot ID        | <a href="#">Q9Y264</a>               |
| Expiration Date   | Please enquire.                      |
| Application Notes | Detection Wavelength: 450 nm         |