

## Product Datasheet

### Human TNNT2 ELISA Kit orb407293

|                   |                                                    |
|-------------------|----------------------------------------------------|
| Description       | Human cardiac isoform of Troponin T,cTnT ELISA Kit |
| Reactivity        | Human                                              |
| Range             | 15.6 pg/ml-1000 pg/ml                              |
| Sensitivity       | 3.9 pg/ml                                          |
| Target            | TNNT2                                              |
| Note              | For research use only.                             |
| Sample Types      | serum, plasma, tissue homogenates                  |
| Uniprot ID        | <a href="#">P45379</a>                             |
| Expiration Date   | Please enquire.                                    |
| Application Notes | Detection Wavelength: 450 nm                       |