

Product Datasheet

Human TNFRSF13B ELISA Kit orb405483

Description	Human Tumor necrosis factor receptor superfamily member 13B(TNFRSF13B) ELISA kit
Reactivity	Human
Range	78 pg/ml - 5000 pg/ml
Sensitivity	19.5 pg/ml
Target	TNFRSF13B
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O14836
Expiration Date	Please enquire.
Application Notes	Detection Wavelength: 450 nm