

Product Datasheet

Human CYB5R3 ELISA Kit orb404893

Description	Human NADH-cytochrome b5 reductase 3(CYB5R3) ELISA kit
Reactivity	Human
Range	25 pg/mL-1600 pg/mL
Sensitivity	6.25 pg/mL
Target	CYB5R3
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P00387
Expiration Date	Please enquire.
Application Notes	Detection Wavelength: 450 nm