

Product Datasheet

Deamidated gliadin peptide antigen orb1925045

Description	This recombinant protein contains multiple deamidated peptide of gliadin with a calculated MW of 37.2 kDa. It is recognized by positive samples in ELISA tests.
Tested Applications	CLIA, ELISA, LF, WB
Form/Appearance	Purified, Liquid
Preservatives	0.1% ProClin300. Tris-buffered saline, pH 8.0
Storage	Upon receipt, quick spin before open vial, aliquot and store at -20°C . Avoid multiple freeze/thaw cycles
Note	For research use only.
Purity	> 90% pure by SDS-PAGE
Expiration Date	6 months from date of receipt.