

Product Datasheet

NPFFR2 Antibody (N-term) Blocking peptide orb1939653

Description	Synthetic peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only.
Protein Sequence	The synthetic peptide sequence used to generate the antibody AP13867a was selected from the N-term region of NPFFR2. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Uniprot ID	Q9Y5X5
NCBI	10886
Expiration Date	6 months from date of receipt.