

Product Datasheet

Scicons Double-stranded RNA (dsRNA) ELISA kit (J2 based) orb1972247

Description	Scicons Double-stranded RNA (dsRNA) ELISA kit (J2 based)
Species/Host	Mouse
Tested Applications	ELISA
Storage	Storage: Upon receipt, store entire kit at -20°C. Once the kit is thawed, you may keep it at 4°C for 5 days. For long-term storage, it is recommended to aliquot and freeze the antibody and dsRNA components at -20°C
Note	For research use only.
Isotype	IgG2a kappa/IgM kappa
Clonality	Monoclonal
Source	The two dsRNA antibodies used in this kit were produced as follows: Female DBA/2 mice were injected intraperitoneally with a mixture of 50 ug L-dsRNA and 75 ug methylated bovine serum albumin, emulsified in complete Freund's adjuvant. After several boosts spleen cells were fused with Sp2/0-Ag14 myeloma cells to generate the hybridoma clones J2 and K2.
Clone ID	J2/K2
Expiration Date	Please enquire.