

Product Datasheet

Cas12 Reaction Buffer, 10X orb1973097

Description	A 10X reaction buffer designed for Cas12 enzyme reactions.
Form/Appearance	Frozen solution of buffer containing: 100mM Tris-HCl, pH 8.0, 500mM NaCl, 100mM MgCl ₂ , 1mg/mL Prionex (0.1% w/v).
Preservatives	Frozen solution of buffer containing: 100mM Tris-HCl, pH 8.0, 500mM NaCl, 100mM MgCl ₂ , 1mg/mL Prionex (0.1% w/v).
Storage	Store at –20°C. To avoid repeated handling and multiple freeze/thaw cycles aliquot product into smaller quantities.
Note	For research use only.
Expiration Date	12 months from date of receipt.
Application Notes	Buffer passed through a 0.2µm filter membrane, packaged aseptically, and confirmed RNase free.