

Product Datasheet

Cas13 Reaction Buffer II, 10X orb1973095

Description	A 10X reaction buffer designed for Cas13 enzyme reactions.
Form/Appearance	Frozen solution of buffer containing: 200mM HEPES, pH 7.0, 500mM KCl, 50mM MgCl ₂ , 1mg/mL Prionex (0.1% w/v)
Preservatives	Frozen solution of buffer containing: 200mM HEPES, pH 7.0, 500mM KCl, 50mM MgCl ₂ , 1mg/mL Prionex (0.1% w/v)
Storage	Store at -20°C. To avoid repeated handling and multiple freeze/thaw cycles aliquot product into smaller quantities.
Note	For research use only.
Expiration Date	12 months from date of receipt.
Application Notes	Buffer passed through a 0.2 µm filter membrane, packaged aseptically, and confirmed RNase free.