

## Product Datasheet

### Recombinant human C22orf32 protein, N-Trx-His orb1974288

Description	Recombinant human C22orf32 protein, N-Trx-His
Endotoxins	Not analyzed
Target	SMDT1
Form/Appearance	Liquid
Concentration	> 0.5mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	N-Trx-His
Note	For research use only.
Protein Sequence	1-64/107
Purity	> 90% as determined by SDS-PAGE
MW	24.5 kDa
Uniprot ID	<a href="#">Q9H4I9</a>
Expiration Date	6 months from date of receipt.
Application Notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.