

## Product Datasheet

### Recombinant human Coagulation factor II protein/F2, C-His (HEK293) orb1974286

<b>Description</b>	Recombinant human Coagulation factor II protein/F2, C-His (HEK293)
<b>Endotoxins</b>	Not analyzed
<b>Target</b>	F2
<b>Form/Appearance</b>	Lyophilized or liquid
<b>Concentration</b>	> 0.5mg/ml
<b>Storage</b>	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
<b>Tag</b>	C-His
<b>Note</b>	For research use only.
<b>Protein Sequence</b>	44-413/622
<b>Purity</b>	> 90% as determined by SDS-PAGE
<b>MW</b>	31.6 kDa
<b>Uniprot ID</b>	<a href="#">P00734</a>
<b>Expiration Date</b>	6 months from date of receipt.
<b>Application Notes</b>	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.