

Product Datasheet

Recombinant HRSVA Major surface glycoprotein G protein (HEK293) orb1974285

Description	Recombinant HRSVA Major surface glycoprotein G protein (HEK293)
Endotoxins	Not analyzed
Target	RSV
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Purity	> 90% as determined by SDS-PAGE
Uniprot ID	P03423
Expiration Date	6 months from date of receipt.
Application Notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.