

Product Datasheet

Recombinant human CD33 protein, C-His (HEK293) orb1974284

Description	Recombinant human CD33 protein, C-His (HEK293)
Endotoxins	Not analyzed
Target	CD33
Form/Appearance	Lyophilized or liquid
Concentration	> 0.52mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only.
Protein Sequence	18-259/364
Purity	> 90% as determined by SDS-PAGE
MW	29.3 kDa
Uniprot ID	Q63994
Expiration Date	6 months from date of receipt.
Application Notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.