

Product Datasheet

SFTPA1 antibody (Biotin) [Out of stock] orb560388

Description	Rabbit polyclonal antivosy against SFTPA1 conjugated to Biotin
Conjugation	Biotin
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Form/Appearence	Lyophilized or Liquid
Concentration	- 100 ?g (in 100 ?l): 1 mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only.
MW	25kDa
Expiration Date	12 months from date of receipt.