

Product Datasheet

Mouse IRAK3(Interleukin 1 Receptor Associated Kinase 3) ELISA Kit orb1994855

Description Mouse IRAK3(Interleukin 1 Receptor Associated Kinase 3) ELISA Kit

Reactivity Mouse

Range 0.32-20 ng/mL

Sensitivity 0.11 ng/mL

Note For research use only.

Sample Types serum, plasma, tissue homogenates and other biological fluids

Expiration Date Please enquire.

Application Notes The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Mouse IRAK3. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Mouse IRAK3. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Mouse IRAK3, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm \pm 10nm. The concentration of Mouse IRAK3 in the samples is then determined by comparing the OD of the samples to the standard curve.