

Product Datasheet

CORONAVIRUS antibody orb433317

Description	Mouse monoclonal antibody to CORONAVIRUS
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	Coronavirus cocktail.
Target	CORONAVIRUS
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Preservatives	
Storage	+4°C, -20°C if preferred
Note	For research use only.
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Clone ID	FIPV3-70
Expiration Date	12 months from date of receipt.